



1. *Chlorophyll a* and *Chlorophyll b* were determined using a spectrophotometer (Shimadzu UV-1601) at 663 nm and 646 nm, respectively. The concentration of chlorophyll was calculated using the following formula: $\text{Chlorophyll concentration (mg/L)} = \frac{\text{Absorbance} \times 1000}{\text{Path length (cm)}} \times \text{Extinction coefficient}$.

NADA QUE MANIFESTAR

722

IMTRA. ADRIANA NOVOA LEDESMA
DIRECTORA GENERAL

J.C. VALERIA JANET HURTADO HURTADO
DIRECTORA ADMINISTRATIVA